

# Fall Prevention Monitors



**Deluxe Bed Sentry**

613-127-1 **\$86.95**



- Dual button reset; must be pressed simultaneously to silence alarm
- Nurse call with remote reset alarm option
- DC adapter jack • No on/off switch



**Deluxe Chair Sentry**

613-129-1 **\$59.95**



- Dual button reset; must be pressed simultaneously to silence alarm
- Nurse call with remote reset alarm option
- DC adapter jack • No on/off switch



**Deluxe Voice Sentry**

613-133-1 **\$65.95**



- Dual button reset; must be pressed simultaneously to silence alarm
- Nurse call with remote reset alarm option
- DC adapter jack • No on/off switch



**Deluxe Pull String**

613-131-1 **\$20.95**

- Pull-string magnet monitor that alarms
- Alarm sounds when the magnetic seal is broken
- Can be integrated into nurse call system
- No on/off switch



**Infrared Bedside Monitor**

613-140-0

**\$51.95**

- Alternative means to detecting residents departure from their bed.
- Emits either a beep or chime alarm when notifying staff
- Multiple installation methods
- Nurse call capability
- Powered by 9volt battery (not included) or optional A/C adapter
- Mounting Hardware included



**Standard Pads**

- Extra large size on bed pads helps to limit false alarms

45-day bed pad	613-128-1	<b>\$23.95</b>
180-day bed pad	613-128-2	<b>\$33.95</b>
1-year bed pad	613-128-4	<b>\$43.95</b>
45-day chair pad	613-130-1	<b>\$17.95</b>
180-day chair pad	613-130-2	<b>\$31.95</b>
1-year chair pad	613-130-4	<b>\$39.95</b>

**Call for complete product line**