

Fall Prevention Monitors



Deluxe Bed Sentry

613-127-1 **\$86.95**

- Dual button reset; must be pressed simultaneously to silence alarm
- Nurse call with remote reset alarm option
- DC adapter jack
- No on/off switch



Deluxe Chair Sentry

613-129-1 **\$59.95**

- Dual button reset; must be pressed simultaneously to silence alarm
- Nurse call with remote reset alarm option
- DC adapter jack
- No on/off switch



Deluxe Voice Sentry

613-133-1 **\$65.95**

- Dual button reset; must be pressed simultaneously to silence alarm
- Nurse call with remote reset alarm option
- DC adapter jack
- No on/off switch



Deluxe Pull String

613-131-1 **\$20.95**

- Pull-string magnet monitor that alarms
- Alarm sounds when the magnetic seal is broken
- Can be integrated into nurse call system
- No on/off switch



Infrared Bedside Monitor

613-140-0

\$51.95

- Alternative means to detecting residents departure from their bed.
- Emits either a beep or chime alarm when notifying staff
- Multiple installation methods
- Nurse call capability
- Powered by 9volt battery (not included) or optional A/C adaptor
- Mounting Hardware included



Standard Pads

- Extra large size on bed pads helps to limit false alarms

45-day bed pad	613-128-1	\$23.95
180-day bed pad	613-128-2	\$33.95
1-year bed pad	613-128-4	\$43.95
45-day chair pad	613-130-1	\$17.95
180-day chair pad	613-130-2	\$31.95
1-year chair pad	613-130-4	\$39.95

Call for complete product line